PROXICLEAN Close proximity cleaning-inplace unit using lost or recovered solutions

AN FEQUANS COMPANY







PROXICLEAN CLOSE PROXIMITY CLEANING-IN-PLACE UNIT



Control hygiene to secure your processes

- > The design of the Proxiclean cleaning station includes all EHEDG, 3A and GMP recommendations
- > The complete automation of these units guarantees the safety, reproducibility and traceability of the cycles as well as the control of all the cleaning parameters.
- Moreover, aware of the environmental impact of cleaning, PIERRE GUERIN provides its expertise to enable sustainable development of your industrial tool by minimizing energy consumption.

Configuration of Proxiclean functions and tanks

It is possible to select the desired functions among the 6 possible (cold rinse, hot rinse, detergent 1, detergent 2, disinfectant, recovery) and to assign each one to a dedicated tank or to use a common tank, according to the following criteria:

Use of one product per dedicated tank

- > The solution is recovered: energy and product savings
- > Shorter cleaning time
- > Higher equipment cost
- Larger overall dimensions

Use of one product per CIP tank

- > The solution is lost Longer
- > cleaning time Optimized
- > equipment cost Smaller
- > overall dimensions

	Dedicated tank	CIP tank 185L on base	In line	From network
Cold water	×			×
Hot water	 Image: A second s	 Image: A second s		×
Detergent 1	 Image: A second s	×	×	
Detergent 2	 Image: A second s	×	×	
Disinfectant	×	×	×	
Recovered water	×	×		

Selection guide

*Several choices are possible

CIP tank		Dedicated tank			
Code	Options	Selection	Code	Options	Selection
R	Clean water tank (CIP tank)	Х	R	Clean water tank (CIP tank)	Х
R1*	Hot water (CIP tank)		Ec*	Dedicated hot water tank	
R2*	Detergent 1 (with dosage) (CIP tank)		D1*	Dedicated detergent tank (with dosage)	
R3*	Detergent 2 (with dosage) (CIP tank)		D2*	Dedicated detergent tank ((with dosage)	
R4*	Recovery (CIP tank)		Re*	Dedicated recovery tank	
R5*	Disinfectant (with dosage) (CIP tank)		Dc*	Dedicated Disinfectant tan (with dosage)	

** only one choice possible

Common selections regardless of the choice of CIP or Dedicated Tank

Code	Options	Selection
P1**	185L vertical tank	
P2**	400L vertical tank	
P3**	800L vertical tank	
lsc1	tank insulation	
Ch1**	Heating with steam	
Ch2**	Heating with hot water	
Ch3**	Heating with electric resistance	
Rc1	1 st rinse with hot water	
Ed**	Direct in-line water	
Ve1**	Water valve on each tank	
Ri1	With feedback sensors switches	
Rp1	With product recovery	
Di1	With in-line disinfection	
Lic1	With continuous level on each tank	
Fte1	With flow meter	
Ftr1	With return flow meter	
FsI1	With return flow sensor	
Tu1	With turbidimeter	
Fta1	With water flow meter	
Sa1	With air blowing	
Pr1	With return pump	
Pe1	With supply pressure sensor	
Fi1	Filter	
Fi2	Return filter	
Ai1	With stainless steel cabinet	

PROXICLEAN CLOSE PROXIMITY CLEANING-IN-PLACE UNIT



Reduced environmental impact

The choice of cleaning devices is decisive in order to combine performance and low flow, water and energy consumption.

CIP balls normal operating conditions	Action radius (m)	Pressure (bar)	Flow rate (m ³ /h)
Static ball (light mechanical action)	1	1.2	1.8
	3	1.2	4.8
	4	1.2	7.2
Rotating ball (mechanical action by sweeping)	1	2	1
	1.5	1.5	5
	2	3	3
Mechanised ball (mechanical action by impact)	2.5	5	3
	4	5	6
	5	5	9
Wall cleaner	2	3	4.5

- > Washing phases with recirculation are preferred to washing with direct delivery of solutions to the effluent
- Stations with recovered solutions integrate a sorting strategy that limits water losses and detergent dilutions
- Recovery of the last rinsing water used is activated as soon as possible
- Tank insulation contributes to reducing energy consumption
- The automatic system allows cleaning times and parameters to be optimised

Cleaning sequence

Test modules allow you to adapt the cycles to your cleaning protocols according to the specificities of your products and equipment.

Type of CIP cycle	Fluids	Temperature	Time (mins)
Pre-wash equipment	Water	15°/ 20°C	5
Blowing & Purging piping	Air		3(1)
Solution preparation detergent 1	Water + caustic solution	60°/ 85°C	7(2)
Recirculation detergent solution detergent 1	Water + caustic solution	60°/ 85°C	15
Rinsing/blowing after detergent 1	Air / Water		3(1)
Solution preparation detergent 2	Water + Acid	50°/ 55°C	7(2)
Recirculation detergent solution detergent 2	Water + Acid	50°/ 55°C	10
Rinsing/blowing after detergent 2	Air / Water		3(1)
Rinsing CIP station and equipment	Water	10°/ 15°C	5
Final rinse equipment	Process water		10
Blowing & Purging piping	Process air		5
Total			60 to 70

(1) Can be carried out as a background task

(2) Can be done as a background task only on a recovered solution station



PIERRE GUERIN SAS

Head office: 179 Grand'Rue BP 40012 79210 Mauzé-Sur-Le Mignon - France Tel. +33 (0)5 49 04 78 00 contact@pierreguerin.engie.com

PIERRE GUERIN Ibérica S.A.U.

c/ Valle de Tobalina, 6 P.I. Villalonquejar 09001 Burgos - Spain Tel: +34 (0)947.259.100

PIERRE GUERIN Deutschland GmbH Daimlerring 6d D-31135 Hildesheim - Germany Tel: +49 5121 708,160 contact.spain@pierreguerin.engie.com contact.germany@pierreguerin.engie.com

PIERRE GUERIN Ltd

3B Swallowfield Courtyard Wolverhampton Road, Oldbury - West Midlands, B69 2JG -United Kingdom Tel: +44 (0) 1452.725.409 contact.uk@pierreguerin.engie.com

PIERRE GUERIN China

Room GHK, 4/F, Building N°3 lane 1505 #100 Qing Yun Road Zhang Jiang High Tech Park Shanghai (201203) - China Tel: +86 (21) 31.33.87.68

contact.china@pierreguerin.engie.com

DCI-Biolafitte

600 North, 54th Avenue ST Cloud, MN 56303 - USA Tel: +1 (800) 671.-7151 info@dci-bio.com





contact@pierreguerin.engie.com www.pierreguerin.fr

FRANCE - UNITED KINGDOM- SPAIN - GERMANY - USA - CHINA



