NEPTUNE™ LAB CONTROL SUITE

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Logiciel Neptune™ Lab

NEW Version 7.0

The screenshot displays the Neptune Lab Control Suite interface, which includes several key components:

- Controller Panel:** A control panel for 'Agit SP (RPM) 1SP1 Ch. 1' with a speed dial ranging from 0 to 1000 RPM. The current value is 525 RPM, and the output is 0.0%. The control state is 'OFF'.
- Process Diagram:** A central schematic of a 200L EVO bioreactor system. It shows various components like the agitator, water inlets (CITY WATER INLET, CLEAR STEAM), and outlets (CITY WATER OUTLET). The diagram is overlaid with a sequence control interface.
- Sequence Control:** A 'SEQUENCE_A' control window showing a table of sequences. The table lists sequence names, vessel numbers, batch numbers, and run times.
- SET_UP CONTROLLER:** A configuration window for the 'Agit SP (RPM) 1SP1 Ch 01' controller. It allows setting parameters such as Control Type (SUP), SetPoint (524.67), Proportional gain (1.0), Integral gain (0), Derivative gain (0), and Cycle Time (5). It also includes options for DDC Out mode, Digital Out Tag, and Vessel Type (Fermentation).
- Digital Output - Actuator:** A control window for the 'IAG1 AGITATOR Ch 01' with a 'Control' button set to 'ON'.
- Sequences Table:** A table listing available sequences:

File	Vessel #	Batch	Step	Run Time
SEQUENCE_A	1	TEST0001	1	0:0:0
SEQUENCE_B	2	12345670	1	0:0:0
SEQUENCE_C	3	decvid19	1	0:0:0
SEQUENCE_D	4	TF123_12	1	0:0:0

TOGETHER BRING LIFE TO YOUR PROJECTS
ENSEMBLE DONNONS VIE A VOS PROJETS

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The **Neptune™ Lab** Software is connected with an optional PLCs to control all phases and process parameters. It also provides the User with an extensive range of supervision and data-logging functions including: Synoptic, Controller parameter editing, Alarms management, Profiles, Trends, Calculations, Strategies, User access, and "Offline Data", through a user friendly and intuitive interface.

This highly flexible and ergonomic system has been designed to control complex processes with single or twin equipment, preferably in Research & Development environments for the Life sciences Industries.

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
- ✓ User friendly and intuitive interface for an installation with single or twin equipment
- ✓ **User Access module** with 3 different levels: operator, engineer and administrator
- ✓ Automatic calibration for the analog inputs and outputs
- ✓ User configurable controllers with selection of inlets / outlets and algorithm type
- ✓ Advanced **Trend, Profile & Offline Data** modules for graphical display of data (Process Values, Set Points, calculation results, Offline Data, etc.) in the forms of trends or historic and controlling set-points via imported or plotted data
- ✓ **Calculation** module to perform automatic calculations (RQ, OUR...) in real-time. The calculation module allows to create formulas and use them with the controllers.
- ✓ **Strategy** module enabling a process to be defined stepwise for process customization and repeatability. Typically a strategy might be used for feed control, complex gassing addition, set-point ramps or event based actions

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
Example of a multiple units architecture

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Example of the **TRYTON™** Line of autoclavable Bioreactors

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The NEPTUNE™ Lab Software offers a user friendly and intuitive interface. A wide range of functional modules are available to suit the specific needs of the Life Science applications.

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